



EUROPEAN SOUTHERN OBSERVATORY

Scientific Report

No. 6 – July 1987

**An Atlas of the
Thorium-Argon Spectrum for the
ESO Echelle Spectrograph
in the $\lambda\lambda$ 3400–9000 Å Region**

S. D'Odorico, M. Ghigo and D. Ponz

Published by
EUROPEAN SOUTHERN OBSERVATORY
Karl-Schwarzschild-Straße 2, D-8046 Garching bei München
Federal Republic of Germany

© Copyright ESO 1987

ESO Scientific Report No. 6

March 1987

An Atlas of the Thorium-Argon Spectrum
for the ESO Echelle Spectrograph
in the $\lambda\lambda$ 3400 - 9000 Å Region

S. D'Odorico, M. Ghigo, D. Ponz



European Southern Observatory
Karl-Schwarzschild-Straße 2
D-8046 Garching bei München

INDEX

Foreword	5
Part 1: Thorium Atlas with 31.6 echelle ($\lambda\lambda$ 3700 - 9000 Å)	7
Part 2: Thorium Atlas with 52.6 echelle ($\lambda\lambda$ 3400 - 4600 Å)	113
Appendix 1: Saturated lines	149
Appendix 2: Line identifications	151

FOREWORD

This atlas is an extended and upgraded version of the "Atlas of the Thorium-Argon Spectrum for the ESO Echelle Spectrograph", ESO scientific report no. 2 in 1984. As the previous one it complements the Operating Manual of the ESO Echelle Spectrograph, CASPEC. It consists of two parts:

1. The first part is based on the CCD frames used for the previous atlas, which cover the spectral range 3700 - 9000 Å. The instrumental set-up was the following: 31.6 lines/mm echelle grating, 300 lines/mm cross disperser, short camera. The Thorium-Argon comparison lamp is from Juniper Co. England. ESO CCD # 3 was used for these exposures. It is an RCA device, with 320×512 pixels of $30 \mu\text{m}$ square.

The following data are presented :

- a) Three plots showing the correlation of echelle orders with wavelength, with resolution and residuals respectively (Figs. 1, 2, 3) and a table showing the variation of the resolution with wavelength and order (Tab. 1).
 - b) The photographic reproductions of 7 CCD frames (Figs. 4-10), taken with CASPEC at different positions of the cross-disperser to include the comparison spectrum from λ 3700 to λ 9000 Å. The orders numbers are also indicated.
 - c) Plots of the individual echelle orders, where the coordinates are the relative intensity versus pixel number. The wavelength identifications are marked on the plots.
2. The second part is based on two new CCD frames and extends the wavelength identifications in the blue-UV side of the spectrum down to 3400 Å. For these exposures the CASPEC set-up was the following: 52.6 lines/mm echelle grating, 300 lines/mm cross-disperser, short camera. The detector was the UV sensitive ESO CCD # 7 (1). It is a GEC CCD with 385×576 pixels, $22 \mu\text{m}$ in size, which has been coated in the ESO lab to enhance the UV sensitivity (2). The data are presented in the same format as in the first part (Figs. 11-15 and Tab. 2).

Appendix 1 is a list of the saturated lines in the red part of the instrumental range that are not indicated in the plots. Appendix 2 contains the identifications of the lines used in the wavelength calibration. The term "W.M." indicates that the listed wavelength is the intensity weighted mean of two or more blended lines. The listing is also available on tape upon request.

In the plots of the spectral orders, the intensities have not been corrected for the blaze function or for variations in the pixel to pixel sensitivity of the CCD.

The accuracy of the wavelength scale has improved considerably with respect to the first edition of the Atlas. This has been made possible by the use of a new, more accurate source for the lines in the Thorium spectrum, the Atlas by Palmer and Engleman (3), in addition to the references used in the first (4, 5, 6, 7, 8). A draft atlas of line identifications

of the CASPEC spectra with the 52.6 lines/mm echelle prepared at La Silla by A. Gilliotte was also useful. The spectral orders were extracted and wavelength calibrated using the Echelle Data Reduction package inside MIDAS, the ESO Data Reduction Software running on the ESO VAX 8600 computers (9, 10). By using a first set of line identifications, a preliminary wavelength calibration was established, leading to additional identifications. The final list is the result of 3 of these iterations. The criteria for identification were the wavelength coincidence, the relative intensities and the energy levels as indicated in the references quoted above. The actual resolution of the plots is at most 2 pixels. Lines separated by less than 2 pixels have not been measured by the automatic procedure and are therefore not included in the Atlas. Fig. 3 and Fig. 13 show the residuals ($\lambda_m - \lambda_o$) for the lines in the final calibration of typical frames taken with 31.6 and 52.6 echelle gratings respectively. For the calibration of the spectra taken with the 31.6 lines/mm echelle, the RMS of the residuals varies from 0.018 Å in the blue to 0.041 Å in the red; this corresponds to about 1/7 of pixel size in Å. For the 52.6 lines/mm echelle the RMS is about 0.013 Å at 4000 Å corresponding to 1/7 of pixel size in Å.

The Atlas has been conceived mainly as an aid to the users of the ESO Echelle Spectrograph. When comparing the plots with spectra taken at the telescope, one should keep in mind that small shifts in the central wavelengths of the orders can be introduced by a change in the inclination of the echelle or a difference in the mounting of the CCD between different runs.

REFERENCES:

1. D'Odorico, S., 1985, The Messenger, No.42, 13.
2. Deiries, S., 1986, Proceedings of the ESO-OHP workshop on the optimization of the use of CCD detectors in Astronomy, J. P. Baluteau and S. D'Odorico editors, published by ESO.
3. Palmer, B. A., Engleman, R. Jr., 1983, Atlas of the Thorium spectrum, Los Alamos National Laboratory. Edited by Helen Sinoradzky.
4. M.I.T. Wavelength Tables, 1969, THE M.I.T. PRESS. (Argon).
5. Moore, C. E., 1959, A multiplet table of Astrophysical Interest, NSRDS-NBS 36 (Argon).
6. Zalubas, R., 1960, New Description of Thorium Spectra, National Bureau of Standards Monograph 17.
7. Junkes, S. J., Salpeter, E. W., 1964, Spectrum of Thorium from 9400 to 2000 Å, Specola Vaticana Publication.
8. Breckinridge, J. B., et al., 1975, Ap. J. Suppl. **29**, 97, Thorium Comparison Spectrum.
9. Ponz, D., Brinks, E., D'Odorico, S., 1986, Proc. SPIE, vol. 627, 707.
10. MIDAS Manual, 1986, ESO Operating Manual No.1.

PART 1

31.6 lines/mm echelle

$\lambda\lambda$ 3700 - 9000 Å




Fig. 1




Fig. 2






Fig. 3

Order	Resolution (Å/pixel)	Central Wavelength (Å)
150	0.117 ± 0.01	3768
140	0.125	4070
130	0.135	4381
120	0.145	4732
110	0.158	5152
100	0.174	5671
90	0.194	6319
80	0.217	7126
70	0.249	8122

Table 1: Resolution for various orders

$\lambda_c = 400 \text{ NM}$



128 →
132 →
136 →
140 →
144 →
148 →
152 →

Fig. 4

$\lambda_c = 475 \text{ NM}$




Fig. 5

$\lambda_c = 550 \text{ NM}$




Fig. 6

$$\lambda_c = 625 \text{ NM}$$




Fig. 7

$\lambda_c = 700 \text{ NM}$




Fig. 8

$$\lambda_c = 775 \text{ NM}$$




Fig. 9

$\lambda_c = 850 \text{ NM}$









Fig. 10
































































































































































































































































104























PART 2

52.6 lines/mm echelle

$\lambda\lambda$ 3400 - 4600 Å





Fig. 11

Fig. 12






Fig. 13

Order	Resolution (Å/pixel)	Central Wavelength (Å)
100	0.078± 0.01	3437
95	0.082	3620
90	0.087	3822
85	0.092	4049
80	0.098	4305
75	0.105	4594

Table .2: Resolution for various orders

$\lambda_c = 350 \text{ NM}$




Fig. 14

$$\lambda_c = 420 \text{ NM}$$









Fig. 15
















































































APPENDIX 1

Saturated Lines

Order	Wavelength (Å)	Order	Wavelength (Å)
82	6965.43	74	7723.76
81	6965.43	74	7724.21
81	7067.22	72	7948.18
81	7068.73	71	7948.18
80	7068.73	71	8006.16
79	7147.05	71	8014.79
78	7272.94	70	8103.69
77	7383.98	70	8115.31
76	7503.87	69	8264.53
76	7514.65	68	8408.21
75	7635.10	68	8424.65
74	7635.10	67	8424.65
		67	8521.45

Table .3: Saturated lines in the spectra (31.6 l/mm echelle)

APPENDIX 2

Line Identifications

wavelength	identification	wavelength	identification	wavelength	identification
3325.120	TH II	3414.46	A II	3480.052	TH I
3327.193	TH I	3415.884	TH I	3480.506	A II
3330.477	TH I	3417.497	TH I	3482.761	TH II
3333.129	TH I	3421.210	TH I	3485.212	TH II
3334.604	TH II	3422.656	TH I	3486.551	TH I
3337.870	TH II	3423.989	TH I	3489.507	TH I
3345.880	TH	3428.999	TH I	3491.536	A II
3348.768	TH I	3433.998	TH II	3493.518	TH II
3350.924	A II	3434.727	TH I	3495.699	TH I
3351.228	TH II	3435.977	TH II	3496.810	TH I
3358.602	TH II	3436.727	TH I	3498.009	TH II
3360.998	TH I	3437.307	TH I	3498.621	TH I
3365.338	TH I	3439.398	TH I	3499.986	TH II
3366.517	TH II	3442.579	TH I	3501.866	TH I
3367.713	TH I-II W.M.	3445.217	TH II	3502.779	TH II
3371.796	TH II	3445.744	TH II	3503.786	TH I
3372.823	TH	3446.547	TH I	3507.546	TH
3373.492	TH I	3449.125	TH I-II W.M.	3509.779	A II
3374.975	TH I	3449.644	TH II	3511.157	TH I
3376.436	A II	3451.702	TH I	3511.562	TH II
3378.573	TH II	3454.095	A II	3514.388	A II
3380.859	TH I	3457.069	TH I	3518.404	TH I
3386.500	TH II	3461.100	TH I-I W.M.	3518.886	TH I
3387.920	TH I	3462.850	TH II	3520.00	A II
3388.531	TH II	3464.127	A II	3521.270	A II
3389.550	TH I-II W.M.	3465.765	TH II	3521.914	TH II
3392.035	TH II	3466.604	TH I-I W.M.	3521.978	A II
3396.727	TH I	3468.219	TH II	3523.506	TH I
3397.516	TH I	3469.345	TH I	3526.634	TH I
3398.544	TH I	3469.920	TH II	3528.411	TH I
3401.711	TH I	3471.218	TH I	3528.954	TH II
3402.695	TH II	3471.959	TH I	3530.515	TH I
3405.558	TH I	3476.747	A II	3531.450	TH I
3408.750	TH I	3478.232	A II	3535.320	A II
3410.075	TH I	3479.172	TH II	3536.011	TH I
3413.013	TH I	3479.684	TH I	3537.159	TH II

wavelength	identification	wavelength	identification	wavelength	identification
3539.587	TH II	3582.355	A II	3621.118	TH II
3542.498	TH I	3583.102	TH I	3622.138	A II
3543.148	A II	3584.175	TH I	3622.795	TH I
3544.018	TH I	3588.441	A II	3623.772	TH I
3545.596	A II	3589.360	TH II	3624.472	TH I
3545.845	A II	3589.749	TH I	3624.895	TH II
3547.338	TH I	3590.925	TH I	3625.628	TH II
3548.514	A II	3591.452	TH I	3626.939	TH I
3549.596	TH I	3592.779	TH I	3629.850	TH
3550.718	TH I	3593.882	TH II	3632.683	A I
3551.402	TH I	3594.985	TH I	3632.830	TH I
3553.110	TH II	3595.617	TH I	3634.582	TH I
3554.306	A I	3597.495	TH	3635.943	TH
3555.013	TH I	3598.120	TH I	3636.834	TH
3555.704	TH I	3599.724	TH I	3637.031	A II
3556.904	A II	3600.432	TH I	3638.319	TH I
3557.464	TH II	3601.034	TH II	3638.644	TH I
3559.508	A II	3603.204	TH II	3639.447	TH II
3559.995	TH I-II W.M.	3604.682	TH	3639.833	A II
3561.030	A II	3605.193	TH I	3642.249	TH I
3561.781	TH I	3605.879	A II	3643.512	TH I
3562.193	A II	3606.522	A II	3645.706	TH I
3563.375	TH I	3608.378	TH	3647.645	TH II
3564.333	A II	3609.445	TH II	3648.421	TH II
3565.030	A II	3610.522	TH I-II W.M.	3649.249	TH II
3567.263	TH I	3612.427	TH I	3649.735	TH I
3569.820	TH I	3612.866	TH I	3649.832	A I
3570.357	TH I	3614.272	TH I-I W.M.	3650.766	TH II
3572.392	TH II	3615.132	TH II	3650.890	A II
3573.219	TH II	3615.850	TH I	3652.168	TH II
3575.322	TH II	3617.117	TH II	3654.461	TH I
3576.557	TH I	3617.671	TH I	3655.278	A II
3576.616	A II	3618.363	TH I	3656.050	A II
3579.336	TH II	3619.212	TH I	3656.693	TH I
3579.457	TH	3620.324	TH II	3657.642	TH I
3581.608	A II	3620.838	TH I	3658.063	TH II

wavelength	identification	wavelength	identification	wavelength	identification
3658.808	TH I	3695.289	TH I	3741.183	TH II
3659.629	TH I	3695.974	TH II	3742.923	TH I
3660.437	TH II	3697.743	TH I	3745.659	TH
3661.621	TH I	3698.106	TH I	3745.970	TH I
3662.750	TH I	3699.181	TH	3747.539	TH II
3663.202	TH I	3699.880	TH I	3749.084	TH I
3663.702	TH II	3700.978	TH I	3749.618	TH I
3666.981	TH I	3703.230	TH I	3751.022	TH I
3668.139	TH	3703.774	TH I	3752.569	TH II
3669.640	TH	3704.861	TH I	3754.593	TH II
3669.968	TH I	3706.767	TH I	3755.212	TH I
3671.539	TH I	3709.862	TH II	3756.294	TH I
3672.300	TH I	3711.304	TH II	3757.694	TH I
3673.26	A	3711.623	TH I	3758.467	TH I
3673.793	TH II	3712.559	TH I	3759.316	TH I
3674.014	TH I	3715.560	TH	3761.469	TH I
3674.903	TH I	3716.583	TH	3762.933	TH I
3675.137	TH I	3718.207	A II	3765.27	A
3675.567	TH II	3719.434	TH I	3766.261	TH
3678.048	TH II	3719.836	TH I	3767.901	TH
3678.27	A	3720.306	TH II	3769.585	TH
3678.480	TH I	3721.215	TH I	3770.056	TH I
3679.134	TH	3721.825	TH II	3770.54	A
3679.710	TH II	3723.656	TH II	3771.370	TH I
3680.061	A II	3724.517	A II	3772.650	TH I
3680.447	TH I	3725.393	TH	3776.271	TH I
3682.486	TH I	3726.724	TH II	3780.966	TH I
3682.545	A II	3727.902	TH I	3783.012	TH II
3684.933	TH I	3729.309	A II	3783.296	TH II
3687.984	TH I	3730.368	TH I	3784.575	TH
3688.760	TH II	3730.748	TH II	3785.600	TH II
3690.623	TH I	3732.985	TH I	3789.167	TH I
3691.411	TH I	3733.672	TH I	3790.355	TH
3691.876	TH I	3734.597	TH II	3790.794	TH I
3692.566	TH I	3737.512	TH I	3792.374	TH I
3694.178	TH	3737.889	A II	3794.698	TH I

wavelength	identification	wavelength	identification	wavelength	identification
3795.385	TH II	3849.183	TH	3915.848	TH I
3797.207	TH I	3850.57	A	3916.417	TH I
3798.103	TH I	3852.135	TH I	3917.269	TH I
3800.198	TH I	3854.510	TH II	3918.070	TH
3801.491	TH I-I W.M.	3856.354	TH I	3919.023	TH I
3803.075	TH I	3859.839	TH II	3920.442	TH
3803.984	TH I	3862.421	TH I	3923.799	TH I
3805.820	TH II	3863.406	TH II	3925.093	TH I
3807.874	TH II	3866.909	TH I	3925.72	A II
3809.49	A	3868.53	A	3927.421	TH II
3810.995	TH I	3869.663	TH I	3928.62	A
3813.067	TH II	3872.722	TH II	3929.669	TH II
3815.566	TH I	3873.148	TH I	3931.256	TH II
3818.685	TH I	3873.474	TH I	3932.55	A
3820.792	TH I	3873.822	TH I	3932.911	TH I
3821.431	TH II	3874.243	TH I	3933.661	TH II
3822.146	TH II	3874.862	TH I	3934.274	TH
3823.067	TH I	3875.373	TH I	3936.353	TH
3825.111	TH I-II W.M.	3877.463	TH I	3937.040	TH II
3825.70	A III	3878.662	TH I	3937.923	TH II
3826.369	TH I	3879.644	TH I	3938.614	TH I
3826.83	A	3884.822	TH	3942.073	TH I
3828.384	TH I	3886.916	TH I	3942.556	TH
3830.060	TH I	3891.97	A	3943.394	TH II
3830.773	TH I	3893.652	TH I	3944.27	A
3831.640	TH I	3895.419	TH I	3945.135	TH
3833.086	TH	3898.437	TH I	3946.10	A II
3834.679	A I	3900.878	TH II	3947.331	TH I
3836.572	TH I-I W.M.	3901.662	TH I	3947.504	A I
3837.875	TH I	3903.102	TH I	3948.030	TH I
3839.700	TH I-I W.M.	3905.186	TH II	3948.979	A I
3840.800	TH I	3908.749	TH I	3950.395	TH I
3841.960	TH II	3910.977	TH II	3951.515	TH II
3842.896	TH I	3911.909	TH I	3952.74	A
3845.42	A	3913.081	TH I	3955.170	TH I
3846.887	TH I	3914.76	A	3955.890	TH I

wavelength	identification	wavelength	identification	wavelength	identification
3956.690	TH II	4009.057	TH I	4063.407	TH I
3959.300	TH I	4009.819	TH I	4064.331	TH I
3960.269	TH I	4011.739	TH I	4067.450	TH I
3961.521	TH	4012.495	TH I	4069.201	TH II
3962.419	TH I	4013.87	A	4069.461	TH I
3964.030	TH I	4014.716	TH I	4070.238	TH I
3967.392	TH I	4017.063	TH I	4072.01	A II
3968.467	TH	4018.099	TH I	4072.40	A
3969.002	TH II	4019.129	TH II	4073.856	TH I
3969.664	TH I	4020.354	TH I	4075.503	TH I
3972.154	TH	4022.076	TH I-II W.M.	4076.64	A
3972.639	TH I	4023.338	TH I	4076.96	A
3973.196	TH I	4024.802	TH I	4079.60	A II
3974.48	A	4025.655	TH II	4080.67	A
3976.414	TH II	4026.147	TH II-II W.M.	4081.367	TH I
3979.36	A	4027.009	TH I	4082.40	A
3980.089	TH	4029.657	TH I	4083.468	TH I
3981.106	TH II	4029.825	TH I	4085.042	TH II
3984.879	TH I	4030.292	TH	4085.434	TH I
3987.206	TH I	4030.842	TH I	4086.520	TH II
3988.027	TH	4032.595	TH I	4087.284	TH I
3988.599	TH II	4033.862	TH I-I W.M.	4088.726	TH I
3990.492	TH I	4035.47	A II	4089.137	TH I
3991.730	TH I	4036.048	TH I	4094.747	TH II
3994.549	TH II	4036.565	TH II	4096.076	TH I
3996.061	TH II	4039.864	TH I	4097.15	A II
3996.668	TH I	4041.203	TH II	4097.747	TH I
3997.865	TH II	4042.91	A II	4100.341	TH I
3998.733	TH I	4043.394	TH I	4102.617	TH I
4000.281	TH I	4044.418	A I	4103.91	A
4001.058	TH I	4048.287	TH I	4105.337	TH II-II W.M.
4001.893	TH I	4049.944	TH I	4107.050	TH
4003.308	TH II	4050.887	TH	4107.853	TH I-II W.M.
4005.092	TH I	4052.94	A II	4108.419	TH II
4007.018	TH II	4053.527	TH I	4109.323	TH I
4008.210	TH I	4059.253	TH I	4110.826	TH I

wavelength	identification	wavelength	identification	wavelength	identification
4112.754	TH I	4181.883	A I	4241.094	TH I
4115.759	TH I	4184.137	TH I	4247.988	TH II
4116.713	TH II	4189.561	TH I	4248.390	TH I
4123.572	TH I-II W.M.	4190.712	A I	4250.314	TH I
4127.412	TH	4191.028	A I	4251.185	A I
4128.65	A	4192.362	TH I	4253.538	TH I
4131.002	TH I	4193.016	TH I	4255.237	TH
4131.712	TH I	4194.936	TH I	4255.751	TH I
4131.73	A II	4198.317	A I	4255.797	TH I
4132.753	TH II	4200.675	A I	4256.253	TH I
4134.068	TH I	4201.846	TH II	4257.496	TH I
4135.480	TH I	4204.041	TH I	4258.520	TH I
4136.285	TH I	4208.411	TH I	4259.361	A I
4136.395	TH	4208.890	TH II	4260.333	TH
4136.436	TH	4210.004	TH I-I W.M.	4266.286	A I
4138.040	TH I	4210.923	TH I	4269.942	TH
4140.235	TH II	4213.067	TH I	4272.168	A I
4142.701	TH II	4214.828	TH I	4272.874	TH I
4143.649	TH I	4216.069	TH I	4273.357	TH II
4148.181	TH II	4217.45	A	4274.024	TH II
4149.987	TH II	4218.69	A	4275.19	A
4154.720	TH I	4220.068	TH I-I W.M.	4276.807	TH II
4156.11	A	4222.67	A	4277.314	TH II
4156.516	TH II	4223.440	TH I	4277.55	A II
4157.280	TH	4223.564	TH	4278.323	TH I
4157.395	TH I	4223.593	TH I	4280.568	TH I
4158.590	A I	4226.299	TH I	4281.067	TH II
4161.738	TH I	4226.726	TH II	4281.414	TH II
4162.509	TH I	4227.02	A II	4282.041	TH II
4164.179	A I	4227.387	TH I	4282.90	A
4165.766	TH I	4228.18	A	4283.009	TH
4168.634	TH II	4229.147	TH I	4283.518	TH II
4170.533	TH I	4229.89	A	4286.228	TH I
4178.059	TH II	4230.426	TH I	4288.668	TH I
4179.253	TH I-I W.M.	4235.463	TH I	4291.809	TH I
4179.715	TH II	4237.23	A II	4294.719	TH I

wavelength	identification	wavelength	identification	wavelength	identification
4297.306	TH I	4348.11	A	4408.883	TH I
4297.99	A II	4349.072	TH I	4412.739	TH II
4299.635	TH I	4352.23	A	4412.893	TH
4299.839	TH I	4352.612	TH I	4414.486	TH I
4300.100	A I	4353.448	TH I	4416.237	TH II
4300.66	A III	4354.482	TH I	4416.844	TH I
4304.956	TH I	4358.320	TH I	4420.90	A
4306.366	TH I	4359.372	TH I	4422.048	TH I
4307.176	TH I	4362.07	A II	4423.720	TH I
4308.122	TH I	4363.794	A I	4426.01	A
4308.600	TH I	4365.930	TH I	4430.18	A
4309.25	A III	4367.87	A II	4431.02	A
4309.940	I-II W.M.	4369.494	TH I	4432.252	TH I
4311.799	TH I	4369.875	TH I	4432.963	TH II
4312.997	TH	4370.76	A II	4433.83	A II
4314.319	TH I	4371.36	A I	4438.746	TH I
4315.254	TH I	4374.124	TH I	4439.123	TH II
4318.415	TH II	4374.785	TH II	4439.48	A II
4320.126	TH II	4375.96	A	4440.324	TH I-I W.M.
4325.274	TH I	4378.176	TH I	4440.866	TH II
4328.915	TH I	4379.74	A	4441.608	TH I
4330.843	TH I	4381.402	TH II	4443.665	TH I
4331.25	A	4381.860	TH II	4445.033	TH I
4332.06	A	4384.656	TH I	4445.315	TH I
4333.560	A I	4385.08	A I	4445.901	TH I
4335.337	A I	4387.787	TH I-I W.M.	4447.834	TH II
4335.709	TH III	4391.110	TH II	4448.554	TH II
4337.277	TH I	4392.974	TH I	4448.88	A II
4338.107	TH I	4393.759	TH I	4450.793	TH I
4340.895	TH I	4397.915	TH II	4452.565	TH I
4342.444	TH I	4400.09	A	4458.001	TH I
4343.951	TH II	4401.02	A	4460.53	A
4345.851	TH I	4401.581	TH I	4461.241	TH I
4344.326	TH II	4402.927	TH I	4461.527	TH I
4345.167	A I	4404.91	A	4463.666	TH I
4346.436	TH I	4408.482	TH I	4465.340	TH II

wavelength	identification	wavelength	identification	wavelength	identification
4469.525	TH I	4533.077	TH I	4611.859	TH II
4470.990	TH I	4533.237	TH I	4612.543	TH I
4474.77	A II	4533.304	TH II	4613.604	TH I
4475.221	TH I	4534.120	TH II	4615.024	TH I
4478.595	TH I	4535.254	TH I	4615.334	TH I
4479.637	TH I	4537.144	TH I	4619.479	TH II
4480.822	TH	4537.67	A II	4621.163	TH I
4481.83	A II	4540.999	TH I	4628.298	TH I
4482.169	TH I	4544.514	TH II	4628.441	A I
4483.346	TH I	4545.08	A	4631.761	TH II
4485.713	TH I	4545.915	TH I	4633.765	TH I
4486.897	TH I	4547.249	TH I	4637.25	A II
4487.495	TH II	4547.78	A	4638.685	TH I
4488.312	TH I	4552.153	TH I	4640.046	TH II
4488.680	TH II	4555.812	TH I	4641.253	TH I
4489.664	TH I	4558.346	TH I	4644.707	TH I
4490.99	A II	4561.347	TH I	4646.686	TH I
4493.333	TH I	4563.661	TH I	4647.251	TH I
4497.914	TH I	4567.240	TH I	4650.234	TH I
4498.55	A II	4568.142	TH I	4651.555	TH II
4498.940	TH I	4570.972	TH I	4651.989	TH II
4499.983	TH I	4577.823	TH I	4655.212	TH
4502.95	A	4579.39	A	4657.94	A
4505.216	TH I	4579.827	TH I	4659.570	TH I
4506.494	TH I-I W.M.	4581.173	TH	4663.202	TH I
4510.526	TH II	4588.426	TH I	4666.005	TH II
4510.733	A I	4589.93	A II	4666.798	TH I
4513.223	TH I	4592.666	TH I	4668.171	TH I
4513.680	TH I	4593.643	TH I	4669.984	TH I
4515.118	TH I	4595.420	TH I	4673.660	TH I
4516.064	TH	4596.097	A I	4676.055	TH I
4519.259	TH I	4598.77	A II	4680.237	TH I
4521.194	TH I	4603.144	TH I	4680.646	TH II
4522.323	A I	4607.934	TH I	4683.351	TH I
4530.57	A II	4608.620	TH I	4686.194	TH I
4532.257	TH II	4609.60	A II	4689.251	TH I

wavelength	identification	wavelength	identification	wavelength	identification
4690.622	TH I	4764.346	TH I	4831.597	TH I
4691.635	TH I	4764.89	A	4832.802	TH II
4694.091	TH II	4765.595	TH I	4833.178	TH I
4695.038	TH I	4766.600	TH I	4840.843	TH I
4695.454	TH I	4773.241	TH I	4843.941	TH I
4700.771	TH I	4774.259	TH I-II W.M.	4845.162	TH I
4702.316	A I	4775.794	TH I	4847.90	A
4703.360	A II	4777.191	TH I	4848.362	TH I
4703.989	TH I	4778.294	TH I	4849.140	TH
4705.760	TH II	4779.728	TH I	4850.439	TH II
4706.251	TH II	4782.761	TH II	4852.868	TH
4708.294	TH I	4783.861	TH I	4858.332	TH II
4712.481	TH I	4786.531	TH I	4861.216	TH I
4712.840	TH I	4787.148	TH I	4861.717	TH I
4720.458	TH	4789.386	TH I	4863.172	TH II
4721.276	TH I	4793.244	TH I	4865.477	TH I
4721.620	A	4795.913	TH I	4865.91	A I
4722.088	TH I	4806.07	A	4867.59	A II
4723.438	TH I	4808.133	TH I	4868.881	TH I
4723.784	TH II	4809.614	TH I	4871.289	TH I
4724.772	TH II	4812.375	TH I	4872.917	TH I
4726.91	A	4813.007	TH I	4874.364	TH I
4728.133	TH I	4813.720	TH	4876.260	A
4729.128	TH I	4813.896	TH I	4878.009	TH I
4732.08	A II	4817.020	TH I	4878.733	TH I
4735.93	A	4818.647	TH II	4879.90	A
4739.676	TH I	4819.193	TH I	4881.204	TH I
4740.529	TH II	4820.464	TH I	4882.25	A
4740.952	TH I-I W.M.	4820.884	TH I	4887.947	A I
4742.117	TH I	4821.587	TH I	4889.06	A
4743.687	TH II	4822.854	TH I	4889.490	TH I
4745.337	TH I	4823.606	TH I	4892.760	TH I
4749.200	TH I	4823.996	TH I	4893.445	TH I
4749.971	TH I	4826.700	TH I	4894.955	TH I
4752.414	TH II	4829.797	TH I	4898.459	TH II
4758.128	TH I-I W.M.	4831.121	TH I	4899.240	TH

wavelength	identification	wavelength	identification	wavelength	identification
4902.054	TH I	4989.308	TH I	5069.338	TH I
4902.77	A II	4993.749	TH I	5081.44	A I
4904.75	A II	5002.097	TH I	5085.02	A II
4910.156	TH I	5003.598	TH I	5085.295	TH
4910.793	TH I	5004.127	TH I	5090.55	A II
4911.378	TH I	5009.35	A	5094.106	TH II
4912.529	TH II	5015.889	TH II	5095.064	TH
4919.815	TH II	5017.16	A II	5096.484	TH I
4920.623	TH I	5017.254	TH II	5098.043	TH II
4921.613	TH II	5019.324	TH II	5098.931	TH I
4925.950	TH	5019.806	TH I	5098.97	A
4927.780	TH I	5022.005	TH I	5100.621	TH I
4929.086	TH I	5028.655	TH I	5101.129	TH I
4933.24	A	5029.891	TH I	5111.278	TH I
4933.852	TH II	5039.230	TH I	5113.382	TH II
4936.774	TH I	5041.122	TH I	5115.044	TH I
4937.829	TH I	5044.719	TH I	5118.200	A I
4939.642	TH I	5045.248	TH I	5122.499	TH I
4943.064	TH I	5047.043	TH I	5125.489	TH I
4945.458	TH I	5048.813	A I	5125.950	TH I
4946.663	TH II	5048.936	TH I	5128.489	TH I
4947.575	TH II	5049.796	TH II	5134.746	TH I
4950.251	TH I	5050.784	TH I	5137.473	TH I
4950.626	TH II	5051.888	TH I	5140.773	TH I
4954.660	TH II	5055.347	TH II	5141.81	A I-II
4960.423	TH	5057.986	TH II	5143.267	TH II
4961.726	TH I	5058.361	TH I	5143.916	TH I
4965.12	A	5059.861	TH I	5145.36	A
4965.731	TH I	5061.656	TH I	5148.211	TH II
4968.755	TH II	5062.07	A	5151.612	TH I
4970.063	TH I-II W.M.	5062.932	TH I	5154.243	TH I
4972.177	TH II	5064.602	TH I	5158.604	TH I
4980.186	TH I	5064.945	TH I	5159.620	TH I
4982.487	TH I	5066.135	TH I	5160.710	TH I-II W.M.
4985.372	TH I	5066.773	TH I	5161.539	TH I
4987.147	TH II	5067.973	TH I	5162.284	A I

wavelength	identification	wavelength	identification	wavelength	identification
5163.458	TH I	5231.159	TH I	5320.770	TH II
5165.82	A I	5233.225	TH II	5322.898	TH I
5168.586	TH I	5234.107	TH I	5325.143	TH II
5168.922	TH I	5237.903	TH II	5326.277	TH I
5173.671	TH I	5238.813	TH I	5326.975	TH I
5175.324	TH I	5239.552	TH I	5329.374	TH II
5175.911	TH I	5240.196	TH II	5330.080	TH I
5176.28	A II	5247.196	TH I	5340.498	TH I
5176.403	TH I	5247.654	TH II	5343.581	TH I
5176.961	TH I	5252.786	A I	5347.971	TH I
5177.623	TH I	5254.259	TH I	5351.126	TH I
5183.989	TH II	5258.360	TH I	5351.836	TH I
5184.453	TH I	5260.104	TH I	5355.636	TH I
5186.413	TH I	5261.472	TH I	5358.706	TH I
5187.337	TH I	5266.710	TH I	5359.827	TH I
5187.467	TH II	5269.792	TH I	5360.150	TH I
5187.746	A I	5273.131	TH I	5361.155	TH I
5190.872	TH II	5274.118	TH I	5362.575	TH I
5193.825	TH II	5277.500	TH II	5369.282	TH I
5194.02	A I	5281.068	TH I	5369.447	TH I
5194.457	TH I	5291.816	TH I	5370.709	TH I
5195.813	TH I	5294.397	TH I	5372.702	TH I
5198.800	TH I	5296.278	TH I	5374.822	TH
5199.163	TH I	5297.743	TH I	5375.352	TH II
5203.848	TH	5298.282	TH I	5376.130	TH
5205.152	TH I	5300.523	TH I	5376.780	TH I
5206.495	TH II	5301.404	TH II	5379.110	TH I
5209.724	TH I	5303.483	TH I	5382.927	TH II
5211.230	TH I	5305.77	A II	5384.301	TH I
5213.349	TH I	5306.986	TH I	5386.610	TH I
5216.596	TH II	5307.465	TH II	5388.050	TH I
5218.527	TH II	5310.266	TH II	5390.440	TH I-II W.M.
5219.109	TH I	5312.002	TH I	5392.572	TH II
5220.926	TH I	5312.528	TH I	5393.972	TH I
5221.270	A I	5312.904	TH I	5394.760	TH I
5228.224	TH I	5317.494	TH I	5398.701	TH I

wavelength	identification	wavelength	identification	wavelength	identification
5398.922	TH I	5488.628	TH II	5571.191	TH I
5399.01	A I	5492.643	TH I	5572.465	TH I
5399.174	TH	5493.204	TH I	5573.353	TH I
5400.145	TH I	5494.330	TH I	5576.204	TH I
5402.69	A	5495.872	A I	5577.70	A III
5403.199	TH I	5499.255	TH I	5579.358	TH I
5407.653	TH I	5499.647	TH I	5580.077	TH I
5410.768	TH I	5501.281	TH I	5580.754	TH I
5415.491	TH I-II W.M.	5504.301	TH I	5583.762	TH II
5417.485	TH I	5506.112	A I	5587.026	TH I
5421.346	A I	5507.538	TH I	5588.750	TH
5421.836	TH II	5508.558	TH I	5593.613	TH
5424.008	TH I	5509.993	TH I	5595.063	TH I
5425.678	TH II	5514.873	TH I	5595.846	TH I
5426.407	TH I	5518.989	TH I	5597.46	A I
5431.112	TH I	5524.582	TH I	5599.654	TH I
5434.151	TH I	5524.93	A I	5601.603	TH I
5435.892	TH II	5527.295	TH I	5602.852	TH W.M.
5437.387	TH II	5528.227	TH I	5604.515	TH II
5439.97	A I	5537.130	TH I	5606.386	TH I
5440.601	TH I	5537.556	TH I	5606.732	A I
5443.118	TH II	5538.608	TH I	5609.573	TH I
5443.21	A I	5539.262	TH I	5610.231	TH I
5447.153	TH	5539.910	TH II	5610.104	TH I
5449.478	TH II	5541.581	TH I	5610.680	TH I
5451.650	A I	5542.890	TH I	5612.068	TH I
5452.218	TH I	5548.175	TH I	5615.319	TH I
5457.37	A I	5551.372	TH II	5619.975	TH I
5458.967	TH	5552.623	TH I	5630.297	TH I
5461.735	TH II	5554.07	A	5632.872	TH II
5464.205	TH I	5557.045	TH I	5633.295	TH I
5470.759	TH I	5558.342	TH I	5639.746	TH II
5473.44	A I	5558.702	A I	5641.734	TH I
5474.865	TH II	5559.891	TH I	5645.668	TH I
5479.075	TH I	5564.201	TH II	5646.451	TH I
5484.146	TH II	5567.998	TH I	5648.991	TH I

wavelength	identification	wavelength	identification	wavelength	identification
5650.703	A I	5771.760	TH I	5871.182	TH I
5657.925	TH I	5772.116	A II	5882.625	A I
5659.130	A I	5773.946	TH I	5885.701	TH I
5664.621	TH I	5777.400	TH I	5886.531	TH I
5665.180	TH I	5781.658	TH I	5888.592	A I
5665.628	TH II	5782.289	TH I-I W.M.	5889.953	TH
5667.128	TH I	5789.645	TH I	5891.451	TH I
5674.986	TH I	5792.430	TH I	5894.698	TH I
5677.053	TH I	5796.068	TH I	5895.281	TH I
5679.005	TH I	5796.418	TH II	5899.844	TH I
5681.900	A I	5797.319	TH I	5904.159	TH II
5685.192	TH I	5798.478	TH I	5905.570	TH I
5691.71	A II	5800.829	TH I	5908.935	TH I-I W.M.
5700.458	TH II	5802.082	A I	5911.229	TH I
5700.917	TH II	5804.141	TH I	5912.084	A I
5702.651	TH I	5807.681	TH I	5914.671	TH
5707.103	TH II	5812.972	TH I	5916.728	TH I
5717.171	TH I	5815.422	TH II	5918.944	TH I
5719.623	TH I	5822.793	TH I	5922.802	TH II
5720.183	TH I	5830.827	TH I	5925.403	TH I
5724.253	TH I-I W.M.	5832.370	TH I	5926.232	TH I
5724.463	TH I	5834.263	A I	5927.13	A I
5725.388	TH I	5838.950	TH II	5928.805	A I
5732.975	TH II	5840.640	TH I	5929.480	TH II
5736.029	TH I	5843.807	TH I	5929.934	TH I
5739.517	A I	5845.918	TH I	5934.409	TH I
5741.170	TH II	5852.680	TH I	5934.461	TH I
5741.829	TH I	5853.474	TH I	5936.386	TH I
5748.741	TH I	5854.120	TH I	5937.162	TH
5749.388	TH II	5859.668	TH II	5937.663	TH I
5753.026	TH I	5860.315	A I	5938.825	TH I
5760.550	TH I	5863.718	TH I	5942.668	A I
5762.794	TH I	5866.811	TH I	5944.647	TH I
5763.529	TH I	5868.374	TH I	5948.799	TH I
5767.778	TH I	5869.850	TH I	5949.26	A I
5768.181	TH I	5870.552	TH II	5955.561	TH I

wavelength	identification	wavelength	identification	wavelength	identification
5959.673	TH I	6049.051	TH I	6125.739	TH II
5962.056	TH I	6052.721	A I	6127.38	A I
5969.737	TH I	6053.380	TH I	6129.545	TH I
5971.59	A I	6055.593	TH I	6137.926	TH I
5973.664	TH I	6059.373	A I	6138.67	A II
5975.064	TH I	6061.536	TH I	6145.43	A I
5981.90	A I	6065.779	TH I	6150.683	TH I
5986.266	TH	6069.020	TH I	6151.993	TH I
5987.289	A I	6073.103	TH II	6154.068	TH I
5989.044	TH II	6077.105	TH I	6154.516	TH I
5991.007	TH I	6077.872	TH I	6155.581	TH I
5994.128	TH I	6078.421	TH I	6157.087	TH I
5996.629	TH I	6079.222	TH I	6161.353	TH I
5999.00	A I	6085.374	TH I	6162.170	TH
6000.762	TH I	6087.262	TH II	6164.479	TH I
6001.203	TH I	6088.030	TH I	6169.822	TH I
6005.165	TH I	6090.76	A I	6172.28	A
6005.74	A I	6098.120	TH	6173.106	A I
6007.072	TH I	6099.083	TH II	6178.431	TH I
6010.160	TH I	6101.16	A I	6180.705	TH I
6011.533	TH I	6101.725	TH I	6182.621	TH I
6013.68	A I	6102.594	TH I	6184.778	TH I-II W.M.
6015.422	TH II	6103.56	A II	6188.125	TH I
6021.035	TH I	6103.641	TH II	6191.905	TH I
6023.224	TH I	6104.568	TH I-II W.M.	6193.856	TH I
6025.14	A I	6105.645	A I	6198.222	TH I
6030.445	TH I	6107.533	TH	6203.492	TH I
6032.124	A I	6112.837	TH II	6205.860	TH I
6032.872	TH I	6114.92	A II	6207.220	TH I
6033.413	TH I	6116.166	TH I	6208.686	TH I
6035.192	TH I	6119.699	TH I	6212.507	A I
6037.697	TH I	6120.556	TH II	6215.945	A I
6038.680	TH I	6121.407	TH	6220.011	TH I
6042.589	TH I	6122.214	TH	6221.319	TH I
6043.230	A I	6123.38	A II	6224.527	TH I
6044.432	TH II	6124.480	TH I	6226.369	TH I

wavelength	identification	wavelength	identification	wavelength	identification
6232.974	TH II	6348.737	TH I	6513.84	A I
6234.855	TH I	6355.910	TH I	6522.043	TH I
6240.953	TH I	6369.139	TH I	6531.341	TH I
6243.13	A II	6369.577	A I	6537.614	TH I
6248.40	A I	6371.943	TH I	6538.115	A I
6250.485	TH I	6376.930	TH I	6545.718	TH I
6257.423	TH I	6379.673	TH I	6551.705	TH I
6258.606	TH II	6384.719	A I	6554.160	TH I
6261.064	TH II	6387.395	TH I	6558.875	TH I
6261.418	TH I	6399.23	A II	6569.632	TH III
6266.173	TH II	6406.446	TH I	6576.122	TH I
6271.544	TH I	6411.899	TH I	6577.214	TH I
6274.116	TH II	6413.614	TH I	6577.656	TH III
6274.164	TH I	6416.315	A I	6580.229	TH I
6277.238	TH II	6431.57	A I	6583.906	TH I
6279.166	TH II	6437.761	TH I	6584.613	TH I
6285.278	TH II	6439.071	TH I	6588.539	TH I
6287.255	TH I	6441.95	A	6591.484	TH I
6291.191	TH I	6443.89	A	6593.462	TH I
6293.242	TH I	6446.771	TH I	6593.939	TH I
6296.876	A I	6450.006	TH I	6596.10	A I
6300.916	TH I	6450.955	TH I	6599.482	TH
6303.250	TH I	6457.282	TH I	6604.853	A I
6307.662	A I	6462.595	TH I-II W.M.	6605.416	TH III
6309.14	A I	6466.56	A I	6613.374	TH I
6310.810	TH I	6483.10	A II	6617.058	TH III
6315.775	TH I	6487.441	TH I	6617.515	TH I
6317.182	TH I	6488.883	TH I	6618.166	TH I
6321.820	TH I	6490.737	TH I	6619.946	TH III
6324.45	A	6493.197	TH I	6638.24	A II
6326.367	TH I	6493.777	TH	6638.912	TH I
6327.277	TH I	6501.992	TH I	6639.72	A II
6331.413	TH I	6503.511	TH II	6643.79	A II
6337.620	TH I	6506.986	TH I	6644.663	TH
6339.668	TH I	6509.050	TH I	6646.540	TH II
6342.859	TH I	6512.364	TH I	6648.495	TH II

wavelength	identification	wavelength	identification	wavelength	identification
6648.958	TH II	6772.174	TH I	6942.538	TH I
6658.677	TH I	6778.312	TH I	6943.610	TH I
6660.64	A I	6780.125	TH I	6945.490	TH I
6662.268	TH I	6780.413	TH I	6948.205	TH I
6664.02	A I	6787.736	TH I	6951.46	A I
6666.36	A II	6788.840	TH I	6954.656	TH I
6668.816	TH I	6791.235	TH I	6965.430	A I
6673.579	TH I	6809.100	TH I	6989.655	TH I
6674.697	TH I	6809.509	TH I	6992.17	A I
6677.282	A I	6818.26	A I	6993.037	TH II
6678.706	TH I	6823.508	TH I	7000.803	TH I
6683.367	TH I	6824.677	TH I	7002.883	TH I
6684.36	A II	6827.240	A I	7018.567	TH I
6687.520	TH I	6829.035	TH I	7030.262	A I
6692.726	TH I	6834.924	TH I	7036.283	TH I
6694.496	TH I	6851.86	A I	7045.797	TH II
6696.140	TH I	6853.523	TH I	7053.619	TH II
6697.712	TH I	6854.109	TH I	7058.489	TH I
6698.85	A I	6861.30	A II	7060.041	TH I
6704.051	TH	6863.52	A II	7060.653	TH I
6711.252	TH I	6866.366	TH I	7064.451	TH I
6713.970	TH I	6866.763	TH I	7067.217	A I
6717.385	TH I	6868.450	TH I	7068.73	A I
6719.20	A I	6871.290	A I	7072.394	TH I
6727.458	TH I	6874.753	TH I	7075.333	TH II
6728.118	TH I	6879.59	A I	7084.169	TH I
6728.759	TH I	6882.811	TH I	7086.70	A I
6733.748	TH I	6886.408	TH I	7089.339	TH II
6742.884	TH I	6886.57	A II	7100.514	TH II
6752.832	A I	6887.10	A I	7107.496	A I
6754.30	A I	6888.17	A I	7109.860	TH I
6756.453	TH I	6889.303	TH II	7114.398	TH I
6757.109	TH I	6909.849	TH II	7122.044	TH
6758.203	TH I	6911.226	TH I	7124.560	TH I
6765.669	TH I	6916.128	TH II	7125.80	A I
6766.56	A I	6937.666	A I	7130.184	TH I

wavelength	identification	wavelength	identification	wavelength	identification
7132.10	TH I	7298.143	TH I	7484.24	A I
7142.331	TH I	7305.404	TH II	7487.973	TH I
7147.041	A I	7308.641	TH I	7489.614	TH I
7148.559	TH I	7311.71	A I	7503.867	A I
7150.284	TH I	7315.066	TH I	7510.42	A I
7154.762	TH I	7316.00	A I	7511.349	TH II
7154.953	TH I	7324.807	TH I	7511.790	TH I
7156.941	TH I	7328.285	TH I	7514.651	A I
7158.83	A I	7329.491	TH I	7523.134	TH I
7159.947	TH I	7335.577	TH II	7525.508	TH II
7162.57	A I	7341.151	TH I	7549.313	TH I
7168.895	TH I	7342.576	TH I	7565.851	TH I
7173.372	TH I	7350.78	A I	7567.741	TH I
7191.133	TH II	7353.316	A I	7569.511	TH I
7200.045	TH I	7372.118	A I	7585.743	TH I-I W.M.
7206.483	TH I	7376.877	TH I	7589.33	A II
7206.986	A I	7383.980	A I	7598.205	TH I
7208.006	TH I	7392.97	A I	7607.823	TH I
7212.689	TH I	7402.252	TH I	7618.33	A I
7218.054	TH I	7411.736	TH I	7625.705	TH I
7219.151	TH I	7412.31	A I	7627.175	TH I
7229.93	A I	7417.790	TH I	7628.86	A I
7230.862	TH I	7418.550	TH I	7630.310	TH I
7233.58	A II	7421.827	TH I	7635.105	A I
7242.092	TH I	7422.26	A I	7647.379	TH I
7244.696	TH I	7425.24	A I	7652.320	TH I
7246.127	TH	7428.940	TH I	7653.828	TH I
7255.354	TH I	7430.253	TH I	7654.699	TH I
7256.986	TH I	7435.33	A I	7658.324	TH I
7258.177	TH I	7436.25	A I	7660.890	TH I
7270.66	A I	7447.801	TH I-I W.M.	7670.04	A I
7272.936	A I	7455.208	TH I	7678.126	TH I
7284.903	TH I	7461.874	TH	7685.307	TH II
7285.44	A I	7471.18	A I	7697.924	TH I
7288.98	TH	7481.354	TH I	7704.81	A I
7296.265	TH I	7483.625	TH I	7710.269	TH I

wavelength	identification	wavelength	identification	wavelength	identification
7713.937	TH I	8079.68	A I	8471.826	TH I
7724.206	A I	8085.219	TH I	8478.358	TH I
7731.738	TH I	8093.623	TH I	8490.30	A I
7782.316	TH	8103.692	A I	8510.624	TH I
7788.934	TH I	8115.311	A I	8516.554	TH I
7798.358	TH I	8122.723	TH I	8521.441	A I
7814.33	A I	8129.405	TH I	8554.944	TH I
7817.770	TH I	8138.475	TH I	8573.120	TH I
7834.457	TH II	8143.138	TH I	8605.78	A I
7840.457	TH I	8159.727	TH I	8620.47	A I
7841.791	TH I	8169.786	TH I	8639.441	TH I
7847.539	TH I	8178.840	A	8667.943	TH I
7861.91	A I	8186.911	TH I	8678.43	A I
7865.969	TH I	8202.147	TH	8709.234	TH I
7868.20	A I	8203.201	TH II	8732.424	TH I
7886.283	TH I	8231.406	TH I	8748.031	TH I
7891.075	A I	8252.393	TH I	8749.169	TH
7900.320	TH I	8254.742	TH I	8758.243	TH I
7916.45	A I	8259.511	TH I	8760.449	TH I
7937.733	TH I	8263.924	TH I	8761.72	A I
7941.726	TH I	8264.521	A I	8766.745	TH I
7948.175	A I	8275.626	TH I	8771.88	A II
7978.973	TH I	8320.855	TH I	8772.805	TH I
7987.973	TH I	8330.450	TH I	8775.573	TH I
8006.156	A I	8332.21	A I	8784.59	A I
8014.786	A I	8358.726	TH I	8841.183	TH I
8022.201	TH I	8367.393	TH I	8849.97	A I
8022.330	TH I	8369.340	TH I	8852.791	TH I
8032.431	TH I	8408.208	A I	8868.833	TH I
8037.23	A I	8416.727	TH I	8892.986	TH I
8046.13	A I	8417.998	TH I	8949.122	TH I
8053.307	A I	8421.225	TH I	8962.19	A I
8062.630	TH I	8424.647	A I	8967.640	TH I
8066.826	TH I	8445.487	TH I	8997.876	TH I
8068.742	TH I	8446.511	TH I	9016.590	TH I
8075.652	TH I	8464.236	TH I	9017.59	A II

wavelength identification

9031.819	TH I
9037.893	TH I
9045.353	TH I
9048.250	TH I
9053.485	TH II
9063.960	TH I
9075.42	A I
9122.966	A I
9194.68	A I
9224.498	TH I
9291.58	A I
9354.218	A I